

# PERFORMANCE & FEATURES

	WIDEX MOMENT™ 440	WIDEX MOMENT™ 330	WIDEX MOMENT™ 220	WIDEX MOMENT™ 110
<b>AUDIBILITY/LISTENING IN QUIET</b>				
Widex PureSound™ with ZeroDelay™ Technology	●	●	●	●
TruAcoustics™ with Personal Gain Integrator	●	●	●	●
Variable Speed Compression	● (IE)	● (IE)	●	●
Digital Pinna	●	●	—	—
Enhanced Dynamic Range Compression	●	●	●	●
Fitting rationale options	WFR/NAL	WFR/NAL	WFR/NAL	WFR/NAL
<b>LISTENING IN NOISE</b>				
True-Input Technology	●	●	●	●
High Definition Locator	15 ch	12 ch	10 ch	6 ch
Speech Enhancer	● (RT/IE)	● (IE)	—	—
Noise Reduction	●	●	●	●
Soft Level Noise Reduction	●	●	●	●
Directional Focus	Via MOMENT App	Via MOMENT App	Via MOMENT App	Via MOMENT App
Reverse focus program	Via MOMENT App	Via MOMENT App	Via MOMENT App	—
<b>USER COMFORT</b>				
Preference Control	●	●	●	●
Fluid Sound Analyser (sound classes)	11 (IE)	7 (IE)	4	3
Multi-directional Active Feedback Cancelling	IE	IE	IE	●
TruSound Softener	●	●	●	—
TruSound AOC	●	●	●	●
Hammer Detector	●	●	●	●
SMARTWIND MANAGER™/Wind Noise Attenuation	●	—	—	—





	WIDEX MOMENT™ 440	WIDEX MOMENT™ 330	WIDEX MOMENT™ 220	WIDEX MOMENT™ 110
<b>USER CONTROLS</b>				
SoundSense Adapt	●	●	●	—
Listening programs	5	4	3	3
InterEar volume/program shift	●	●	●	●
Speech (S)/Music (M) streaming template	S/M	S/M	S/M	—
Phone+	●	●	—	—
SmartSpeak/SmartTone	●	●	●	●
Personal Zen	●	●	●	●
Zen/Zen+	IE	IE	IE	IE
Wireless CROS compatible	●	●	●	●
<b>HEARING HIGH FREQUENCIES</b>				
High-Frequency Boost	●	—	—	—
Audibility Extender	●	●	●	●
<b>FITTING AND FINE TUNING</b>				
Processing channels	15	12	10	6
Fine-tuning channels	15	12	10	6
High-resolution Sensogram (9 fitting bands)	●	●	●	●
Sound Diary/Datalogging	●	●	●	●
<b>APPS, ACCESSORIES AND REMOTE CONTROLS</b>				
WIDEX MOMENT™ App with SoundSense Learn*	●	●	●	●
WIDEX TONELINK™ App	●	●	●	●
COM-DEX App	●	●	●	●
TV PLAY™**	●	●	●	●
COM-DEX / COM-DEX REMOTE-MIC	●	●	●	●
RC-DEX	●	●	●	●
TV-DEX	●	●	●	●
FM+DEX	●	●	●	●

\* mRIC R D and RIC 312 D only, no SoundSense Learn in 110

\*\* Available with WIDEX MOMENT™ compatibility late 2020.



# WIDEX MOMENT™ COMPARISON CHART

WIDEX MOMENT™ hearing instruments are available in a variety of technology levels, each designed to accommodate your specific listening requirements.

WIDEX MOMENT™ 440	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★
WIDEX MOMENT™ 330	★★★★	★★★★	★★	★★★★	★★★★	★★	★★★★	★★★★	★★★★	★★★★	★★
WIDEX MOMENT™ 220	★★★	★★★	★★	★★★	★★★	★★	★★★	★★★★	★★★	★★★	★★
WIDEX MOMENT™ 110	★★	★★	★★	★★	★★	★★	★★	★★	★★	★★	★★

WIDEX MOMENT™ hearing aids are packed full of the latest technology to provide the most natural sound ever heard in a hearing aid.


	BENEFIT / USE	440	330	220	110
SMARTWIND MANAGER™	Makes speech and other sounds audible in wind - enjoy conversations in the great outdoors.	●	-	-	-
High Frequency Boost	Clear and crisp sound in high frequencies (in music, for example).	●	-	-	-
Digital Pinna	Mimics the normal pinna effect to help locate where a voice is coming from and hear from behind.	●	●	-	-
Speech Enhancer	Reduces background noise and enhances speech - follow the conversation even in noisy environments.	●	●	-	-
Phone+	Cancels the mic on your opposite ear when using the phone to reduce background noise. Great for landline phones or when streaming is unavailable.	●	●	-	-
Soundsense Learn / Adapt	Widex MOMENT™ App uses artificial intelligence to guide you to optimal, personalised hearing. Remembers your preferences for better sound in every listening situation.	●	●	●	-
HD Locator	Better speech understanding in noise.	15ch	12ch	10ch	-
Sound Classes	Automatically adapts to different listening environments, for optimal clarity, comfort and audibility.	11	7	4	3
Programs	Select the programs you would like your hearing aid to have and manually adjust to appropriate program for different environments.	5	4	2	2

# WIDEX MOMENT™ SOUND CLASS TECHNOLOGY QUICK GUIDE

Exclusive Sound Class Technology allows **WIDEX MOMENT™** hearing aids to adapt automatically to different situations - delivering optimal clarity, comfort and audibility in every environment.

SOUND CLASS	LISTENING ENVIRONMENT	USE / BENEFIT	440	330	220	110
<b>MUSIC</b>	Listening to music	Enjoy your favourite music as it's meant to be heard.				
<b>MUSIC:</b> Classical or contemporary	Optimises sound for different types of music	Settings alter for different levels of intensity and variations in sound level / frequency, associated with different music genres.				
<b>QUIET</b>	Quiet environments with limited background noise (e.g. library)	Hear softer background noises for a more natural hearing experience.				
<b>QUIET</b> with Speech	Optimises speech sounds in quiet environments	Hear conversations in quiet environments comfortably, even where speech is soft or hushed.				
<b>PARTY</b>	Where there is lots of background noise / conversation (e.g. Party)	Distracting background noises are less audible.				
<b>PARTY</b> with Speech	Optimises speech sounds in noisy social situations	Focus on the conversation you want to hear in a noisy social situation.				
<b>TRANSPORT</b>	Situations with low-frequency sounds (e.g. car noise)	Reduces audibility of continuous low-frequency sounds.				
<b>TRANSPORT</b> with speech	Optimises speech sounds in situations with low-frequency sounds	Hear the conversation when travelling.				
<b>URBAN</b>	Where there are variable sounds (e.g. supermarket, city, noisy workplace)	Enjoy the sounds of the city naturally and comfortably without becoming disorientated or overwhelmed.				
<b>URBAN</b> with speech	Optimises speech sounds in urban environments	Comfortably hear the conversation without over amplification of variable background noises.				
<b>SOCIAL</b>	Where multiple conversations occur simultaneously (e.g. busy restaurant)	Small social gatherings where multiple people are talking around you. Focus on the conversation, without distraction.				

## WIDEX MOMENT™ KEY FEATURES





THERE'S NEVER BEEN  
A BETTER MOMENT

WIDEX **MOMENT**<sup>TM</sup>

**WIDEX**

SOUND LIKE NO OTHER



## WHAT MAKES A GREAT HEARING AID?

Designed to help you participate more fully

Promotes hearing well to keep your mind sharp

Produces natural sound to make hearing easier

Provides flexible hearing in any environment

Customized solution to set you up for success



# HEARING LIKE YOU REMEMBER

Until now, even the best hearing aids sounded like hearing aids. That's all changed with the enhanced Widex Moment with PureSound™ technology. Unlike other devices, Moment processes sound with lightening speed, so there is none of the delay that causes that annoying “tinny” artificial sound. Just pure, natural sound so you hear every moment, with sound that over 90% of wearers are satisfied with.<sup>1</sup>



“WHEN I AM WEARING  
MY MOMENT HEARING AIDS,  
I FORGET I HAVE THEM IN.”

—BARB G

# WIDEX SOUNDRELAX SOOTHING SOUNDS AIMS TO HELP YOU REDUCE STRESS

In today's busy world, we are all becoming increasingly aware of the effects stress can have on our minds and bodies. This is why finding new and innovative ways to overcome anxiety and to relax has never been more important.

**Did you know our Widex Moment hearing aids can help you regain your inner peace?**

Widex SoundRelax™ provides a library of sounds to choose from to help relaxation, well-being, and concentration. This unique feature is designed to help all hearing aid wears — with or without tinnitus — deal with anxiety, soothe the mind and boost concentration. SoundRelax, when activated in your hearing aid, is accessible in the Widex Moment app.



Download the Widex Moment app  
to get started

\* The Widex MOMENT app works with all 2.4 GHz compatible Widex MOMENT hearing aid models. Use Widex Tonelink app for all other models.





## MORE WAYS TO CONNECT

Widex Moment offers direct streaming from iPhone® and compatible Android™ devices so you can stream high-quality calls, video, and music directly to your hearing aids.

With the Moment app, you have full control of your hearing aids, so it's never been easier to fine-tune your hearing to get exactly the sound you prefer — in the moments that are important to you.

### **BETTER HEARING MADE EASY**

With MyGuide, a feature in the Moment app designed to help you get the most out of your hearing aids, find step-by-step instructions, videos, and short courses that will help you with everything from turning your hearing aids on to personalizing their sound.

### **SOUND MADE FOR YOU**

The MySound™ 2.0 section of the app uses the power of AI to help adjust your sound. You can choose to have the AI guide you through a series of sound comparisons to tailor the sound to your particular situation. Or you can select a sound profile that has been automatically generated from the preferences of other Moment users in similar situations all over the world!



Download the Widex Moment app  
to get started

## BUILT TO LAST

Widex Moment is designed for maximum moisture resistance and long life — because less corrosion means more durability.

The advanced, water-resistant nano coating and the fact that our rechargeable options have no battery door means there is less opportunity for moisture to sneak in — providing you with a listening experience free from disruptions.

## DISCREET RECHARGEABILITY

Convenience is discreet with Widex Moment. Although the sRIC R D model is the market's smallest rechargeable receiver-in-canal hearing aid, you can rest easy knowing you'll have all the power you need to hear every moment.

Take better care of your Widex Moment sRIC R D hearing aids with Widex Charge n Clean. The charger offers:

- four-hour charge for a full day use
- 30-minute charge for four hours of use
- Sanitizing function of earware
- Drying function

Widex Charge n Clean





## WIDEX MOMENT MODELS AND STYLES

Widex Moment is available in a variety of colors to match your skin tone, hair color, or personality. Below are a few of our most popular options. Your Hearing Care Professional will help you find the perfect fit.

### HEARING AID MODELS



sRIC R D

#### RECEIVER-IN-CANAL (RIC)

The case sits behind the ear with the receiver in the ear-tip instead of the case. This allows for a smaller design without sacrificing quality.

- sRIC R D is the smallest rechargeable hearing aid on the market with a modern design
- For mild to severe hearing loss



BTE R D

#### BEHIND-THE-EAR (BTE)

The case sits behind the ear with a hook over the top of the ear. Hearing aids behind the ear are the perfect choice if you want extra power in a stylish design.

- Rechargeable options available
- T-coil comes standard
- For mild to severe-to-profound hearing loss



CIC

#### IN-THE-EAR (ITE)

The custom hearing aid solution sits discreetly in your ear. This solution is available in a variety of models from the easy handling ITE to the ultra small completely-in-the canal hearing aid.

- Customized to fit your ear canal
- For mild to moderately severe hearing loss

# WIRELESS ACCESSORIES

Accessories help you effortlessly connect to audio equipment, remote controls, television and mobile or landline phones. Ask your Hearing Care Professional for more details.

## BLUETOOTH CONNECTIVITY



reddot award 2018  
winner

### SOUND ASSIST\*

Multi-functional device featuring remote mic, hands-free, t-coil and remote control

### TV PLAY™\*\*

Stylish and stable TV streaming

## ADDITIONAL WIRELESS CONNECTIVITY\*\*



### COM-DEX

For hands-free streaming



### COM-DEX REMOTE MIC

For improved audibility in noise



### RC-DEX

Your remote control



### PHONE-DEX 2\*\*\*

For land line phone streaming



### TV-DEX

For enjoying TV




### FM+DEX

For flexible streaming

\* For compatible styles sRIC R D, RIC 312 D, BTE R D & BTE 13 D

\*\* For all styles with WidexLink™ technology (excludes CIC-M).

\*\*\* Additional handset sold separately. Must have a PHONE-DEX 2 to use additional handset.





# WIDEX **MOMENT** APP

QUICK GUIDE

FOR ALL WIDEX **MOMENT**™ BLUETOOTH HEARING AIDS

**WIDEX**

SOUND LIKE NO OTHER



## GETTING STARTED



The Widex Moment app is designed to work with Widex Moment Bluetooth® hearing aids.

With the app, you can control and personalize your hearing aids and stream sound directly from iOS and compatible Android devices.

You can download the Widex Moment app from the App Store or Google Play.





# PAIRING WIDEX **MOMENT** WITH YOUR PHONE

## APPLE iOS

1. Enable Bluetooth by going to Settings > Bluetooth > On
2. Go to Settings > Accessibility > Hearing Devices
3. Place your hearing aids in pairing mode by turning them off and back on again
4. Select your hearing aids from MFi hearing devices and press Pair for each hearing aid
5. Your hearing aids are now paired with your phone and are Moment ready to stream music and calls
6. Open the app to connect your hearing aids to the app

## ANDROID

1. To enable Bluetooth, swipe down from the top of your phone screen and select the Bluetooth icon
2. Open the Moment app and select Get Started
3. Choose the hearing aid model, Rechargeable or with Batteries
4. Place your hearing aids in the pairing mode by turning them off and back on again
5. Then choose Pair for each hearing aid
6. Setup is now successful, press OK

# WIDEX MOMENT APP HOME SCREEN

In-app help

Connection  
or battery status

Reset sound

More options  
(for non-sound  
adjustments)

Save changes as  
personal program


Mute and volume  
General sound  
adjustments

Select program

Home and Menu

My Sound

My Guide




STAY  
UPDATED

You will occasionally receive automatic notifications when updates are available for your phone and your hearing aids.

It's important that you accept the updates to have the most recent functions installed.

You can get information and tutorials about how to get the most out of your MOMENT hearing aids via the My Guide section of the app.



## SHARE YOUR REAL-LIFE DATA

In the Widex Moment app, you'll be asked to give consent to sharing some data.

First, you'll be asked to share data with Widex. This helps us improve your listening experience and make future hearing aids even better.

Next, you may also be asked to give your consent to sharing data with your hearing care professional. You will only be asked if you've already agreed to do this with your hearing care professional.

Sharing these data means your hearing care professional will get access to the hearing adjustments you make with the app. That way, they can give you an even more personalized fitting.

All data are securely encrypted.

## DESIGN YOUR OWN HEARING

You can always adjust the sound in your hearing aids. Your adjustments won't change the settings your hearing care professional has created.

You can adjust the sound in three ways:

- My Sound including SoundSense Learn
- General sound adjustments
- Directional focus

# SOUND MADE FOR YOU

The My Sound section of the app harnesses the power of artificial intelligence (AI) to help you adjust your sound to suit you.

You can choose to have the AI guide you through a series of sound comparisons so you can tailor the sound to your particular situation, or you can select a sound profile that has been automatically generated from the preferences of other Widex Moment users in similar situation all over the world!

1. Enter **My Sound**
2. Press **Get Started** to Optimize your sound
3. Select a **Situation** and **Intention**
4. The app offers you two AI based **Made for You** recommendations...
5. Or choose **Create Your Own** to have our AI guide you to a more customised sound via a series of sound comparisons.

You can access the feature from the bottom of the home screen.

# CREATE YOUR OWN

When you select Create Your Own, you are using SoundSense Learn. SoundSense Learn is a simple task that uses AI to easily guide you to a sound that's personalized for specific situations.

SoundSense Learn presents two sound profiles, A and B. Listen to the two sound profiles and use the blue slider to indicate the one you prefer. When you press **Next**, two new sound profiles that are based on your choice appear.

We recommend that you continue the process until completed. Then you can save your adjustments as a personal program.



# GENERAL SOUND ADJUSTMENTS

Press the  icon on the home screen to access general sound adjustments:

- L /R volume: Adjust the volume for each hearing aid individually
- Equalizer: Adjust the sound in Bass, Middle or Treble to get a more personal sound – or choose one of the pre-sets at the bottom of the screen
- Sound mixer: Choose whether you want your hearing aids to focus on your surroundings or a second input, like streaming (only visible during streaming and specific programs)

# DIRECTIONAL FOCUS

With directional focus you can choose to focus on sound that's coming from the front, back, left or right.

Select the directional focus program  from your programs on the home screen and choose the direction you want to focus on.

# PERSONAL PROGRAMS

When you've made your adjustments and have found a sound that you're satisfied with, you can save it. Choose **Save as** on the home screen and give the program a name and icon. You can also add a location to the program by long pressing the program. That way, the hearing aids automatically switch to that program when you revisit the location.



# TROUBLE SHOOTING

## If you experience problems with pairing:

- Check that Bluetooth® is activated and that the hearing aids are close to your phone
- Change/charge the hearing aid batteries
- Turn your hearing aids off and then on again. The hearing aids are in pairing mode for the first three minutes after they have been turned on
- Deactivate and activate Bluetooth® in your phone's settings
- For iOS, make sure that the Widex Moment app has access to Bluetooth®: Go to Settings > Moment and activate Bluetooth®

## If it still doesn't work:

- Unpair and re-pair your hearing aids
- Close the Widex Moment app if it's active
- Delete the MOMENT app and reinstall it
- Fully power off your phone and turn it back on
- Ensure all phone updates are installed

## If Bluetooth is interfering:

- Move away from any potential Bluetooth interference
- Reduce the distance/barriers from the phone to the hearing aids
- Decrease the number of apps/functions occurring on the phone
- Move your phone out of your pocket or purse

## If you experience streaming problems:

- Make sure the hearing aids are paired with the phone
- For iOS, go to Settings > Accessibility > Hearing Devices > Audio Routing. Make sure that **Call Audio** and **Media Audio** are both set to **Automatic** or **Always Hearing Devices**
- Deactivate and activate Bluetooth in your phone's settings
- Disable Wi-Fi on your phone
- Verify that your phone has all the current updates installed
- Close all open apps that are not in use
- Turn AirDrop or Nearby Share setting to OFF if not in use

## Get more information

<https://www.widex.com/en-us/hearing-aids/apps/moment-app/>


## Is your phone compatible?

<https://www.widex.com/en-us/hearing-aids/apps/moment-app/>

## Want tips or want to learn more about hearing and hearing loss?

<https://www.widex.com/en-us/blog>

# HEARING LOSS AND DEMENTIA



ASK US HOW WIDEX **MOMENT™** CAN HELP.

**WIDEX**

SOUND LIKE NO OTHER



# WIDEX SOUND CLASS TECHNOLOGY

Widex EVOKE hearing aids are available in a variety of technology levels, each designed to accommodate individuals' specific hearing needs. Whether a person's hearing loss is mild, severe, or somewhere in between, there is an EVOKE hearing aid that will meet their needs, budget, and lifestyle.

Our exclusive Sound Class technology helps EVOKE hearing aids adapt automatically to different situations, delivering optimal clarity, comfort, and audibility. EVOKE hearing instruments categorize listening environments into the eleven distinct sound classes illustrated below. Although they adapt automatically to each listening environment, Hearing Healthcare Professionals can further refine these settings to deliver optimal audibility and comfort to every EVOKE wearer.

SOUND CLASSES	SOCIAL	MUSIC	QUIET	QUIET W/ SPEECH	PARTY	PARTY W/ SPEECH	TRANSPORT	TRANSPORT W/ SPEECH	URBAN	URBAN W/ SPEECH
<b>EVOKE 440</b> 11 Available Sound Classes		 CLASSICAL and CONTEMPORARY								
<b>EVOKE 330</b> 7 Available Sound Classes	Unavailable	 MUSIC						Unavailable		Unavailable
<b>EVOKE 220</b> 4 Available Sound Classes	Unavailable	 MUSIC		Unavailable		Unavailable		Unavailable	Unavailable	Unavailable
<b>EVOKE 110</b> 3 Available Sound Classes	Unavailable	 MUSIC		Unavailable		Unavailable	Unavailable	Unavailable	Unavailable	Unavailable

## HOW DO SOUND CLASSES HELP EVOKE WEARERS HEAR BETTER?

Automatic SoundClasses are available in all levels of EVOKE and allow for additional clarity and comfort. Indicates additional speech prioritization.

<b>Social</b>	Hear your best in small groups with multiple speakers: meetings, classes, family dinners
<b>Music Classical and Contemporary</b>	Hear your best with music by emphasizing the differences in intensity, frequency, and variations between classical and contemporary music
<b>Quiet</b>	Hear your best in quieter situations: home, doctor visits, 1:1 meetings, golf, and quiet office settings
<b>Party</b>	Hear your best in noisy situations like restaurants, parties, and large groups
<b>Transport</b>	Hear your best while on the move in cars, buses, and trains
<b>Urban</b>	Hear your best in a variety of environments: city streets, shopping malls, supermarkets, and noisy office settings

Learn how Widex technology features adapt to suit the various sound classes.



## SOUND CLASSIFICATION **DYNAMIC FEATURE ASSESSMENT**

Not all sound classes are available in all hearing aids.

SOUND CLASS	BALANCED / DEFAULT	MORE AUDIBILITY	MORE COMFORT
 <b>Social</b>	<ul style="list-style-type: none"> <li>Optimized speech understanding when there are multiple speakers</li> <li>Focus on speech coming from multiple directions</li> <li>Maximum audibility of speech sounds</li> <li>Attenuates distant soft sounds</li> </ul>	<ul style="list-style-type: none"> <li>Speech intelligibility optimized for improved speech understanding when it is in close proximity to the listener</li> <li>Maximum contrast between loud and soft sounds</li> </ul>	<ul style="list-style-type: none"> <li>Attenuation of overall sound level</li> <li>Speech intelligibility optimized for comfort</li> <li>Maximized impulse sound control</li> </ul>
 <b>Music Classical</b>	<ul style="list-style-type: none"> <li>Preserves full dynamic range of sound</li> <li>Increased audibility of faint passages</li> <li>Sounds heard from all directions</li> <li>Fast compression</li> <li>Speech in noise off</li> <li>Wind noise reduction off</li> <li>No impulse control</li> </ul>	<ul style="list-style-type: none"> <li>Minimal impulse control</li> <li>Attenuation of soft mid-frequency sounds</li> </ul>	<ul style="list-style-type: none"> <li>Attenuation of overall sound level</li> <li>Provides less overall gain</li> <li>Maximized impulse sound control</li> </ul>
 <b>Music Contemporary</b>	<ul style="list-style-type: none"> <li>Preserves full dynamic range of sound</li> <li>Steady gain</li> <li>Maximizes naturalness</li> <li>Full low frequency sound</li> <li>Sounds heard from all directions</li> <li>Slow compression</li> <li>High-frequency boost (6K Hz +)</li> <li>Speech in noise off</li> <li>Wind noise reduction off</li> <li>No impulse control</li> </ul>	<ul style="list-style-type: none"> <li>Minimal impulse control</li> <li>Minimal attenuation of mid-frequency sounds</li> <li>Increase gain for mid- and soft-level low and high frequency sounds</li> <li>Provides more overall gain</li> </ul>	<ul style="list-style-type: none"> <li>Attenuation of overall sound level</li> <li>Provides less overall gain</li> <li>Maximized impulse sound control</li> </ul>
 <b>Quiet</b>	<ul style="list-style-type: none"> <li>Audibility of soft speech</li> <li>Reduction of soft background noise</li> <li>Stable sound picture/less compression</li> <li>Adaptive sound directionality</li> <li>Attenuates impulse sounds to maintain comfort</li> <li>Minimal noise management</li> <li>Wind noise attenuation</li> </ul>	<ul style="list-style-type: none"> <li>Increased audibility for all soft sounds</li> <li>More overall gain and loudness</li> </ul>	<ul style="list-style-type: none"> <li>Maximum reduction of soft level sounds</li> <li>Less overall gain and loudness</li> <li>Maximized impulse sound control</li> </ul>
 <b>Party</b>	<ul style="list-style-type: none"> <li>Optimized speech understanding in noise</li> <li>Focus on speech in front</li> <li>Maximum audibility of speech sounds</li> <li>Maximum attenuation of low frequency noise</li> <li>Attenuates impulse sounds</li> <li>Wind noise attenuation</li> </ul>	<ul style="list-style-type: none"> <li>Speech intelligibility optimized for improved speech understanding</li> <li>Provides more overall gain</li> <li>Increased contrast between loud and soft sounds</li> </ul>	<ul style="list-style-type: none"> <li>Speech intelligibility optimized for comfort</li> <li>Provides less overall gain</li> <li>Maximized impulse sound control</li> <li>Reduced contrast between loud and soft sounds</li> </ul>
 <b>Transport</b>	<ul style="list-style-type: none"> <li>Maximum attenuation of low frequency noise</li> <li>Adaptive sound directionality</li> <li>Attenuates impulse sounds</li> <li>Stable sound picture</li> <li>Active speech and noise attenuation</li> <li>Access to speech signal</li> <li>Wind noise attenuation</li> </ul>	<ul style="list-style-type: none"> <li>Speech intelligibility optimized for improved speech understanding</li> <li>Minimal impulse control</li> <li>Provides more overall gain</li> <li>Increased contrast between loud and soft sounds</li> </ul>	<ul style="list-style-type: none"> <li>Speech intelligibility optimized for comfort</li> <li>Provides less overall gain</li> <li>Maximized impulse sound control</li> <li>Reduced contrast between loud and soft sounds</li> </ul>
 <b>Urban</b>	<ul style="list-style-type: none"> <li>Audibility for soft, moderate and louder sounds</li> <li>Adaptive sound directionality</li> <li>Minimum attenuation of low frequency noise</li> <li>Maintains a wide input dynamic range</li> <li>Wind noise attenuation</li> <li>Attenuates impulse sounds</li> </ul>	<ul style="list-style-type: none"> <li>Speech intelligibility optimized for improved speech understanding</li> <li>Provides more overall gain</li> <li>Minimal impulse control</li> </ul>	<ul style="list-style-type: none"> <li>Maximized impulse sound control</li> <li>Provides less overall gain</li> <li>Reduced contrast between loud and soft sounds</li> </ul>